

REPORT DOCUMENTATION PAGE				Form Approved OMB No. 0704-0188	
Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing this collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to Department of Defense, Washington Headquarters Services, Directorate for Information Operations and Reports (0704-0188), 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to any penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number. <b>PLEASE DO NOT RETURN YOUR FORM TO THE ABOVE ADDRESS.</b>					
1. REPORT DATE (DD-MM-YYYY) 07-06-2010		2. REPORT TYPE Conference Proceeding		3. DATES COVERED (From - To) 2009-2010	
4. TITLE AND SUBTITLE  Enhanced Interferometry with Programmable Spatial Light Modulator				5a. CONTRACT NUMBER NNX08CA25C	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S)  James D. Trolinger and Joshua S. Jo				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES)  MetroLaser, Inc. 8 Chrysler Irvine, CA 92618				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING / MONITORING AGENCY NAME(S) AND ADDRESS(ES)  NASA Goddard Space Flight Center				10. SPONSOR/MONITOR'S ACRONYM(S) NASA GSFC	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION / AVAILABILITY STATEMENT  DISTRIBUTION STATEMENT A: Approved for public release; distribution is unlimited.					
13. SUPPLEMENTARY NOTES Presented at Mirror Technology Days, Boulder, Colorado, USA, 7-9 June 2010.					
14. ABSTRACT Modern optical components such as aspheres are difficult to inspect. Interferometry, the usual solution, lacks the required dynamic range. Too many fringes are produced. Adapting interferometry for null measurements currently necessitates special optical components for each inspected component. In this effort, a programmable Spatial Light Modulator (SLM) was incorporated into an optical inspection instrument to boost the dynamic range and enable null testing on almost any component using off the shelf optics. The instrument was used to measure the surface accuracy of a 6 inch tall cone with a top diameter of 8.4 inches and a bottom diameter of 8.2 inches. Measurement simulations were compared with test data and the system errors were investigated. A procedure was also developed and tested to subtract system aberrations. A dynamic range of 150 wavelengths was demonstrated.					
15. SUBJECT TERMS Aspheres, Interferometry, Spatial Light Modulator, Surface Accuracy, Optics, Mirror, Zernike, Freeform Optics, Null Testing, Hartman, Wavefront					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT  SAR	18. NUMBER OF PAGES  43	19a. NAME OF RESPONSIBLE PERSON Hans-Peter Dumm
a. REPORT UNCLASSIFIED	b. ABSTRACT UNCLASSIFIED	c. THIS PAGE UNCLASSIFIED			19b. TELEPHONE NUMBER (include area code) 505-853-8397

# Enhanced Interferometry with a Programmable Spatial Light Modulator

Presented at  
Mirror Technology SBIR/STTR Workshop  
June 7<sup>th</sup> to 9<sup>th</sup>, 2010  
Millennium Hotel, Boulder CO

by  
James D. Trolinger and Joshua S. Jo  
MetroLaser, Inc.  
jtrolinger@metrolaserinc.com

This work was supported by NASA Goddard Space Flight  
Center under SBIR Phase II contract #NNX08CA25C  
(POC's: Geraldine Wright and David Content)

*This presentation is approved for General Public release*



# Presentation Summary

- Summarize the problem & innovation
- Describe the hardware and software
  - Digital Interferometry with preconditioned wavefronts
  - Hybrid Hartmann/Digital Interferometry
- Typical measurements
- Potential Applications
- Future Work



# The Problem Being Addressed

Modern optical components such as aspheres are difficult to inspect. Interferometry, the usual solution, lacks the required dynamic range. Too many fringes are produced. Adapting interferometry for null measurements currently necessitates special optical components for each inspected component.



# Innovation/Solution


- Incorporate a programmable SLM into an optical inspection instrument enabling:
  - Combined interferometry and Hartmann Sensing
  - Virtually unlimited dynamic range
  - Extended dynamic range of interferometry by preconditioning waves
  - Null testing on almost any component using off the shelf optics.
- Incorporate Instantaneous Digital Interferometry Technology



# **Challenges Faced and Solved**

- **Incorporating the SLM into a PhaseCam**
- **Integrating Hartmann and Digital Interferometry**
- **Calibrating the SLM**
  - **Phase-only mode by controlling polarization**
  - **Corrected gamma curve, i.e. linear phase shift versus grayscale value**
- **Impressing the required phase function on the SLM**
- **Identifying and minimizing errors and noise**

# Conventional Shack-Hartmann Characteristics



- Lenslet diameters,  $d$ , define spatial resolution over the wavefront being measured.
- $\delta$  (sensitivity) proportional to  $f$ , which should be less than  $d$  to prevent confusion

Without aberration



With aberration



'Measuring aberration of the eye with wavefront technology'

Giuseppe Colicchi, et al

Zur Veröffentlichung eingereicht bei "Physics Education", 2006


# Key Components: SLM, Pixelcam\*

## SLM for a compensator

- Programmable holographic optical element
- Produce wavefronts of any shape and can simulate freeform optical surfaces
- Holoeye SLM
  - SLM can produce a phase up to  $2\pi$  at 632.8nm
  - Assign 0 to 255 grayscale values to 0 to  $2\pi$  (or  $1\lambda$ )
  - Can generate higher phases by wrapping phase
  - Can provide more than 150 wave tilt
  - Pixel size  $\sim 19$  microns.
  - 19.5 x 14.6mm size (1024 x 768 pixels)

## Pixelated Phasecam for a detector

- Spatial phase shifting interferometer
- Single shot, insensitive to vibration



0	90
270	180


Phase information of the object,  $\Delta\phi(x, y)$  can be obtained from the 4 intensities on each unit cell.

$$I(x, y) = \frac{1}{2} (I_r + I_s + 2\sqrt{I_r I_s} \cos(\Delta\phi(x, y) + 2\alpha_p))$$


\*Produced and Trademark by 4D Technologies, Inc, Tucson, AZ <http://www.4dtechnology.com>



# Hybrid Hartmann & Digital Interferometer




# Using a SLM as a Scanning Shack-Hartmann Component




# Scanning Shack-Hartmann System

Control Screen




Control screen is a map of a scanned aperture (or pencils light) of the SLM.

Data Screen




**No Overlap of Focused Spot  
Programmable Aperture Size**



# Scanning Shack-Hartmann System



Control Screen



'Pencils' of light are deflected by the test object's slope of wavefront as they scan the object.


**SLM corrects Aberrations**

Data Screen




# Scanning Shack-Hartmann System

Control Screen




Data screen maps, in time, the sequence of angular deviations caused by the object.

Data Screen



**Dynamic range is limited only by the CCD size.**




# Constraints in conventional Hartmann Testing that can be obviated with this concept

- Sensitivity requires a longer focal length which can cause focused image overlap of adjacent light pencils.  
**> SLM enables scanning in time, NO OVERLAP**
- If the wave is aberrated, the focused spot will not be round, so its centroid is more difficult to locate, further reducing accuracy  
**> SLM enables preconditioning/aberration correction**
- Spatial resolution is limited to the diameter of the lenses in the array.  
**> Limited to fractional pixel size of SLM**




# Mandrel Used for Demonstration

Photo of the Mandrel




Geometry of the Mandrel





- Mandrel provided by NASA GSFC, 6-inches tall, top and bottom diameters are different, cone shape.

# System Design (Top view)





# The Hardware





6-7-2010-Monday




# Simulation Results

Beam size at cylinder lens:  $\sim 1.8$  inch  
150mm fl cylinder lens


Top view




Side View



Interferogram




# Reducing the number of Fringes from a mandrel




(Various waveplates and telescopes not shown)

Measured Interferogram



Applied Phase on the SLM



Residual Interferogram









# Compensation of one of Mandrel's wavefronts for a one inch beam

## Before Compensation




## After Compensation




# Resolution & Signal to Noise

-Tilt by the flat mirror,  $\sim 7$  waves




Tilt introduced by flat mirror



Interferogram

$\sim 14$  fringes



3D view

PV.  $\pm 3.4 \lambda$




- Residual wavefronts, less than  $0.05 \lambda$



A 3D visualization of a seismic volume. The volume is shown as a rectangular prism with a top surface displaying seismic data. The x-axis is labeled 'x-axis' and the y-axis is labeled 'y-axis'. The z-axis represents depth, with labels on the left side: 0, 250, 500, 750, 1000, 1250, 1500, 1750, 2000, 2250, 2500, 2750, 3000, 3250, 3500, 3750, 4000, 4250, 4500, 4750, 5000, 5250, 5500, 5750, 6000, 6250, 6500, 6750, 7000, 7250, 7500, 7750, 8000, 8250, 8500, 8750, 9000, 9250, 9500, 9750, 10000. The x-axis has labels 0, 250, 500, 750, 1000, 1250, 1500, 1750, 2000, 2250, 2500, 2750, 3000, 3250, 3500, 3750, 4000, 4250, 4500, 4750, 5000, 5250, 5500, 5750, 6000, 6250, 6500, 6750, 7000, 7250, 7500, 7750, 8000, 8250, 8500, 8750, 9000, 9250, 9500, 9750, 10000. The y-axis has labels 0, 250, 500, 750, 1000, 1250, 1500, 1750, 2000, 2250, 2500, 2750, 3000, 3250, 3500, 3750, 4000, 4250, 4500, 4750, 5000, 5250, 5500, 5750, 6000, 6250, 6500, 6750, 7000, 7250, 7500, 7750, 8000, 8250, 8500, 8750, 9000, 9250, 9500, 9750, 10000. The volume is colored with a grayscale gradient, with darker colors representing higher values. The top surface is labeled 'S010'.



6-7-2010-Monday

# Determining System Aberrations with a Cat's Eye



**Due to SLM, cylinder lens, optical alignments, etc**


Coaxial position






# Cat's eye wavefronts by ZEMAX based on Experimental results

Experimental measurement




Simulated Fringes



Leftover <  $0.1\lambda$



Fitted Zernike Polynomials upto 36 terms, found and put into Zemax




Simulated Cats' eye wavefronts by ZEMAX














## Simulated, No System Aberrations

## Measured



Beam : ~27mm




Wavefronts at cat's  
eye position at each  
height

Wavefronts at  
coaxial position at  
each height


Wavefronts at  
coaxial position at  
each height



# Net fringe due to a Mandrel




**Coaxial position**



**Cat's eye position**



**Simulated Net Fringes**




**Experimental Net Fringes**




# Wavefronts from the Mandrel, II

Simulated Wavefronts  
of the Mandrel



Experimental Wavefronts





# Unique Features & Applications

- Very wide dynamic range
- Applicable to aspheres and non axis-symmetric optics, i.e. freeform optics
- Enables null point testing
- Enables removing system aberrations

# Summary and Conclusions

- Demonstrated wide dynamic range digital interferometry/Hartmann for advanced optical components using an SLM to:
  - provide wavefront preconditioning.
  - Hartmann & Interferometry in the same instrument
  - Scanning Hartmann
  - extended dynamic range.
  - Null point testing.
- Hartmann provides information needed to program SLM for wavefront preconditioning
- SLM Calibration procedure to produce Gamma curve.
- Procedure to subtract system aberrations.
- Accuracy of  $\lambda/20$  with SLM in system.
- Concept extends dynamic range by more than  $150\lambda$




# Future Work

- Software Development
  - Transforming Hartmann data into wavefront Preconditioning data.
  - Automating Calibration
  - Interferogram stitching
  - System automation
- Hardware improvements
  - Calibration to push system accuracy
  - Incorporating improved SLM's
  - System packaging

**Back up Slides**




# Including System Aberrations



Beam : ~27mm


Simulated



cat's eye position

coaxial position

Measured





coaxial position



# Key Components: SLM, Pixelcam\*


## SLM for a compensator

- Programmable holographic optical element
- Produce wavefronts of any shape and can simulate freeform optical surfaces
- Holoeye SLM
  - SLM can produce a phase up to  $2\pi$  at 632.8nm
  - Assign 0 to 255 grayscale values to 0 to  $2\pi$  (or  $1\lambda$ )
  - Can generate higher phases by wrapping phase
  - Can provide more than 150 wave tilt
  - Pixel size  $\sim 19$  microns.
  - 19.5 x 14.6mm size (1024 x 768 pixels)



## Pixelated Phasecam for a detector

- Spatial phase shifting interferometer
- Single shot, insensitive to vibration



0	90
270	180

Phase information of the object,  $\Delta\phi(x, y)$  can be obtained from the 4 intensities on each unit cell.

$$I(x, y) = \frac{1}{2} (I_r + I_s + 2\sqrt{I_r I_s} \cos(\Delta\phi(x, y) + 2\alpha_p))$$

\*Produced and Trademark by 4D Technologies, Inc, Tucson, AZ <http://www.4dtechnology.com>

# Enhanced Interferometry with a Programmable Spatial Light Modulator

*MetroLaser Incorporated  
Irvine, CA*

## **INNOVATION**


A unique optical inspection system incorporates a dynamic holographic optical element that combines and extends both Hartmann and digital interferometry with preconditioned wavefronts. The resulting system exhibits a wide dynamic range and will be especially useful for inspecting aspherical and free form optics

## **ACCOMPLISHMENTS**


- ❖ Demonstrated wide dynamic range optical inspection system
  - Wavefront preconditioning via spatial light modulator.
  - Hartmann & Interferometry in the same instrument
  - Scanning Hartmann
  - Extended dynamic range  $>150\lambda$ .
  - Null point testing.
- 
- ❖ Hartmann provides information for programming SLM/HOE
  - ❖ Procedure to subtract system aberrations.
- 
- ❖ Accuracy of  $\lambda/20$  with SLM in system.

## **COMMERCIALIZATION**

- ◆ Basic concept patent application has been submitted
- ◆ Marketing to manufacturers of free form and aspherical optics
- ◆ Currently beta testing systems in service tests
- ◆ A wide range of customers have indicated interest
- ◆ Alliances made with two specific manufacturers
- ◆ This technology can provide significant savings of time and money over competitors
- ◆ This technology can enable inspections not provided by any competitors



*Before Correction*




*After Correction  
Wavefront Preconditioning*

## **APPLICATIONS/FUTURE WORK**

- ◆ Government and Commercial Applications
  - X-ray telescope mirrors and mandrels
  - Free form Optics
  - Aspherical Optics
  - SiC telescope Mirrors
  - Ogive windows
- ◆ Future Work
  - Additional Software for System Operation and Automation
  - System Packaging
  - Noise floor and resolution improvement

# SLM Calibration

SLM



Gray scale value is varied on upper half of SLM

Gray scale value is held constant at 0 on lower half of SLM

Applied gray values on the SLM

Gray value: 63

...

Applied gray values on the SLM


Gray value: 127

Gray value: 0


Gray value: 0

Interferogram Result

Phase shift vs. gray scale value is measured interferometrically




S: Phase shifted amount  
P: Period of fringes





- Relative phase shift recorded to achieve an accuracy of data
- Gray value was varied from 0 to 255 on the top half, while it was held constant (0) on the bottom half
- Phase shifted vs applied gray value on the SLM produces gamma curve

# Corrected Gamma Curve




The calibrated gamma curve (blue color), is used to linearize the phase vs. gray scale response (red curve)

*Errors when producing  $1\lambda$  of tilt and corrected one after calibration*





$S = 0.25\lambda$  shift when 63 value applied

# Corrected Gamma Curve



The calibrated gamma curve (blue color), is used to linearize the phase vs. gray scale response (red curve)

*Errors when producing  $1\lambda$  of tilt and corrected one after calibration*




$S = 0.25\lambda$  shift when 63 value applied



# Producing cylindrical wavefronts


## Applied Phase on the SLM



$$\Phi \sim a y^2$$

*Matlab formula used to generate the phase.*

## Interferogram on the pixelcam




- Cylindrical wrapped phase written by Matlab code and applied on the SLM, which is good candidate compensating phase for a Mandrel.
- Phase produced and shown by the interferogram on the pixelcam



# Producing arbitrary wavefronts

## Interferogram on the pixelcam


### Applied Phase on the SLM




$$\Phi = 20 \cdot (y - 0.25)^3 - x^2 + x \cdot y + x \cdot y^3$$

Matlab formula used to generate the phase

- Arbitrary wrapped phase written by Matlab code and applied on the SLM
- Phase produced and shown by the interferogram on the pixelcam




### 3D view




# Compensation of Tilt

-Applied compensating tilt by the SLM


Wrapped phase applied on the SLM



(Various waveplates and telescopes not shown)




~ 14 fringes




Compensating wrapped phase applied on the SLM




# Compensation of Cylindrical wavefronts




Simulation




Measured Interferogram on the pixelcam



Applied Phase on the SLM





Residual fringes on the pixelcam




# Compensation of Cylindrical wavefronts


## Extreme Case




Measured Interferogram



Applied Phase on the SLM









Residual Interferogram



The phase on SLM is used to bring the dense fringes down into the measurement range of the pixelcam. i.e: extended dynamic range.



# Wavefront Differences at Different Heights of Mandrel



# Enhanced Interferometry with a Programmable Spatial Light Modulator

*MetroLaser Incorporated  
Irvine, CA*

## **INNOVATION**


A unique optical inspection system incorporates a dynamic holographic optical element that combines and extends both Hartmann and digital interferometry with preconditioned wavefronts. The resulting system exhibits a wide dynamic range and will be especially useful for inspecting aspherical and free form optics

## **ACCOMPLISHMENTS**


- ❖ Demonstrated wide dynamic range optical inspection system
- Wavefront preconditioning via spatial light modulator.
- Hartmann & Interferometry in the same instrument
- Scanning Hartmann
- Extended dynamic range  $>150\lambda$ .
- Null point testing.
- ❖ Hartmann provides information for programming SLM/HOE
- ❖ Procedure to subtract system aberrations.
- ❖ Accuracy of  $\lambda/20$  with SLM in system.

## **COMMERCIALIZATION**

- ◆ Basic concept patent application has been submitted
- ◆ Marketing to manufacturers of free form and aspherical optics
- ◆ Currently beta testing systems in service tests
- ◆ A wide range of customers have indicated interest
- ◆ Alliances made with two specific manufacturers
- ◆ This technology can provide significant savings of time and money over competitors
- ◆ This technology can enable inspections not provided by any competitors



*Before Correction*



*After Correction  
Wavefront Preconditioning*

## **APPLICATIONS/FUTURE WORK**

- ◆ Government and Commercial Applications
  - X-ray telescope mirrors and mandrels
  - Free form Optics
  - Aspherical Optics
  - SiC telescope Mirrors
  - Ogive windows
- ◆ Future Work
  - Additional Software for System Operation and Automation
  - System Packaging
  - Noise floor and resolution improvement